

Technical Datasheet

Analysis Name: Brix Determination of Citrus Juices and Concentrates

Method Number: NQA-00.0925, NQA-00.0926

Scope of Application: This method determines the acid corrected Brix solids of citrus

juices and concentrates having less than or equal to 30.2%

titratable acidity as anhydrous citric acid.

Description: The refractive index (Brix) of a substance is the ratio of the

velocity of light in a vacuum to its velocity in the substance; this, in turn, is dependent on composition, concentration, and temperature of the substance. For practical measurements, the scales of standard refractometers indicate refractive indices

with respect to air rather than vacuum.

The Brix scale corresponds to the change of refractive index of a sucrose solution in water as a function of the sucrose concentration. The 'Brix value is equal to the % (mass/mass) of sucrose. For liquids which are mostly water and sugar, such as non-citrus juices, the °Brix value is a close approximation of the % dissolved solids. For liquids with a high proportion of nonsugar dissolved solids, this approximation is less accurate. Citrus juices are such liquids, having a high proportion of dissolved organic acids, mostly citric acid.

The uncorrected °Brix value is first measured using a refractometer. The titratable acidity value is also measured. A correction factor in °Brix is chosen from a table based on the titratable acidity result and is added to the uncorrected °Brix value. This is the acid corrected Brix value and is a better

approximation of the total dissolved solids.

Sample Weight 25 q **Required:**

Analytical Platform: Automatic Refractometer, pH Meter, Autotitrator

Special Information: Juice

Single Strength calculation included with BRIX_SS variation

6625 Eiterman Rd, Dublin, OH 43017

P (614)526.5200

E ngacdublininfo@us.nestle.com

TDS-NQA-00.0926-2



Analyte Reported	Alias	Unit of Measure	Limit of Quantification	Reproducibility
Brix		ºBx	N/A	°Bx<10, r= 0.2 °Bx ABS; 10≥°Bx<30, r=2% RSD; °Bx≥30, r=1%RSD
Brix, Acid Corrected		ºBx	N/A	N/A
Brix, Single Strength		ºBx	N/A	N/A