



NQAC

Nestlé Quality Assurance Center
Dublin

Technical Datasheet

Analysis Name: Insoluble Solids (Pulp Determination) in Citrus Juices and Concentrates

Method Number: NQA-06.6030

Scope of Application: This method applies to the determination of percent pulp in fruit juice concentrates and purees.

Description: The sample for analysis is centrifuge, and the content of pulp is determined as percentage by volume of the juice.

Sample Weight Required: 50 g

Method Reference: Determination of Centrifugable Pulp, International Federation of Fruit Juice, Procedures #60 (1991)

Analytical Platform: Automatic Refractometer / Centrifuge

Special Information: Juice

Analyte Reported	Alias	Unit of Measure	Limit of Quantification	Reproducibility
Pulp		%	4	25%