



NQAC

Nestlé Quality Assurance Center
Dublin

Technical Datasheet

Analysis Name: L-Malic Acid Enzymatic Determination

Method Number: NQA-06.8113

Scope of Application: Determination of L-Malic acid in juice concentrates, juice purees, juice pulps and other similar products.

Description: Rapid determination of L-Malic acid in juice concentrates and juice purees via enzyme addition and detection by spectrophotometric absorbance.

Sample Weight Required: 50 g

Analytical Platform: UV-Spectrophotometer

Analyte Reported	Alias	Unit of Measure	Limit of Quantification	Reproducibility
Malic Acid		mg/100 g	2	10%