



NQAC

Nestlé Quality Assurance Center  
Dublin

# Technical Datasheet

**Analysis Name:** Coffee Brew Standard

**Method Number:** Coffee Brew Standard

**Scope of Application:** Coffee Brews

**Description:** Coffee Grounds and/or Pods may be submitted for in-house brew prior to analytical testing.

**Sample Weight Minimum:**  
**Required:** 45g per Brew (for Drip Machines)  
# Brews + 1 (for Pods)  
**Recommended:**  
85g per Brew (for Drip Machines)  
# Brews X 2 (for Pods)

**Analytical Platform:** Sample Prep

**Special Information:** Customer must populate special instructions to indicate brewer, any specific instructions for brewing, and number of brews. Standard brewing parameters are stated in the table below.

Brewed coffee may be stored up to 5 days for analysis. Prices may fluctuate based on difficulty.

Brewer	Number of Brews	Volume of Water (or Final Brew Weight)	Weight of Ground Coffee
<b>Drip Machine</b>	1-5	720 mL	40 g
<b>Keurig</b>	1-5	8 oz	N/A
<b>Other</b>	1-5	Specify in Special Instruction	N/A