

## **Technical Datasheet**

Analysis Name: Total, Soluble, & Insoluble Dietary Fiber in Foods

Method Number: OM-AOAC-991.43\_SDF\_IDF

**Scope of Application:** This method applies to processed foods, grain and cereal

products, fruits and vegetables which do not contain low molecular weight fiber, or in regions where low molecular weight fiber is excluded from the definition of total dietary fiber.

Not applicable to oily samples.

**Description:** Duplicate test portions of dried foods undergo sequential

enzymatic digestion by heat stable a-amylase, protease, and amyloglycosidase to remove starch and protein. For total dietary fiber (TDF), enzyme digestate is treated with alcohol to precipitate soluble dietary fiber before filtering, and TDF residue is washed with alcohol and acetone, dried, and weighed. For insoluble and soluble dietary fiber (IDF and SDF), enzyme digestate is filtered, and residue (IDF) is washed with warm water, dried and weighed. For SDF, combined filtrate and washes are precipitated with alcohol, filtered, dried, and weighed. TDF, IDF, and SDF residue values are corrected for

protein, ash, and blank

Sample Weight 25 g

Required:

Method Reference: AOAC 991.43

**Analytical Platform:** Enzymatic, Gravimetric

Special information: If only total dietary fiber results are needed, OM-AOAC-

991.43 TDF should be chosen.

This method is not suitable to determine non-digestible oligosaccharides (NDOs) such as fructoligosaccharides (FOS), inulin, galactoligosaccharides (GOS), resistant maltodextrins, polydextrose. Low molecular weight fiber should be analyzed

with AOAC 2017.16 or AOAC 2011.25.

Presence of gum (ie guar gum) needs to be communicated at

the time of sample submission.



Analyte Reported	Alias	Unit of Measure	Limit of Quantification	Reproducibility
Total Dietary	TDF	g/100 g	0.5	Culinary: 3-60%
Fiber				Cereals: 3-50%
Insoluble Fiber	IDF	g/100 g	0.5	Culinary: 3-60%
				Cereals: 3-50%
Soluble Fiber	SDF	g/100 g	0.5	Culinary: 3-60%
				Cereals: 3-50%